
ACCESS **B4**

& **B4** **XL**

Manual de **ES**
instrucciones

Instruction manual **EN**

ACCESS
by **BLAUTEC**®

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

Description

Fixed hydraulic lift ACCESS B4


The purpose of the ACCESS lifts is to facilitate access to the pool for people with reduced mobility. It requires a minimum depth of 800 mm and supports a maximum weight of 150 kg.

Smooth and progressive functioning

The elevator moves vertically, making a rotation around its own axis of approximately 160°. The force for the movement is extracted from the water network pressure itself. It should be between 3.5 and 5 bars of pressure. The lift is composed of a tubular structure in stainless steel 316L anchored to the pool by three plugs and three adjustable supports.

Technical characteristics

	B4	B4XL
Max. lifting capacity:	150 Kg	
Rotation angle:	160°	
Required pressure:	3,5 a 5 bars	
Anchor arm:	355	655 mm
Minimum depth:	800 mm	
Weight:	43 Kg	44,5 Kg
Dimensions:	2425 x 700 mm	
Structure:	Stainless Steel AISI 316L	
Chair:	Synthetic material	

Pre-installation

It is important to be aware of the following before installing the lift:

- The installation must be carried out by personnel with technical knowledge.
- Installation requires more than one operator (+40 Kg).
- Consider the current accessibility regulations according to your region.
- Choose a good location for the elevator, always maintaining a minimum distance of 1 meter from any architectural element of the pool (stairs, skimmers, lights, etc.) and ensure a minimum depth of 800 mm.

We recommend placing the elevator in shallow areas so that, if necessary, a caregiver or assistant can help the user with greater difficulties from inside the water.


- Have a nearby water outlet that provides the necessary pressure (3.5 and 5 bars). If the pressure is lower, a pump should be installed to increase pressure. If, on the other hand, it is higher, a pressure reducer will be needed.

- Check the **packaging contents**:

- ACCESS lift
- Valve operation shaft ²² and prisoner ⁴⁰
- 3 Metallic nuts D12x100 ¹³
- Cream for steel maintenance
- 2 Allen wrenches (2,5 y 4)
- Instruction manual

Anchor dimensions

The anchor will be located behind the overflow grille. The size of the anchor arm of ACCESS B4 and B4 XL allows to cover most pool profiles.

Assembly instructions

The steps to follow for the installation of the lift are outlined below.
Remember to perform all steps with the lift empty (without load).

What do you need?

Hammer drill &
12 mm drill bit

Hammer

Allen
wrenches
(included)


Fixed
wrenches 13 (2)
and 17

Level


Level

1. Place the lift in the pool. The anchor leg ¹⁰ must be positioned after the overflow grille at an angle between 105° and 110° with respect to the body of the lift.


2. Adjust holding shaft ¹¹ to **ensure the verticality of the elevator** using the level for that purpose. Adjust the position of lower swivel ¹⁹, leaving it fully supported against the pool wall. If necessary, to achieve the level, you can adjust the height or angle of the leg of the anchor support assembly ⁷.

IMPORTANT: do not loosen the set screws on the elements that comprise the unit ^{3 2}, as this will dislocate the entire cylinder ¹ assembly and damage its joints.


3. With the elevator leveled, proceed to **mark the anchoring holes** ¹³ of the anchor plate on the floor.

4. Remove the elevator and proceed to **drill the floor** with a 12 mm drill bit. Ensure the depth of the holes so that the anchor screws enter completely. Clean the dust thoroughly.


5. Place the elevator and the three anchoring plugs ¹³, plus the ground terminal, tightening them until they are securely fastened. Check the correct position of the plugs.

6. Once secured, **recheck the vertical level** and proceed with corrections, if necessary, with the lower joint and the anchor support assembly (steps 2 and 3).


7. **Check the bolts and screws of the elevator** to ensure they are all properly tightened.


8. **Place the chair fixing assembly** ⁹ in position with the M8x16 nuts ³⁵ with washer ³⁶ and their female counterparts ³⁷, to the support flanges assembly ⁸.

9. Position the supplied prisoner ⁴⁰ and shaft ²² to **connect the upper drive pieces** ²¹ with the upper connecting flanges assembly ³.


10. **Connect the lift to the water supply** in order to supply the circuit. The lift has a special connector enabling quick and easy connection ²³.

11. Proceed with its start-up. **Perform 3 or 4 cycles** to purge any air that may remain in the sleeve. Verify that the sleeve rises slowly and without fluctuations.

12. **Check the position of the seat.** If necessary, its height or rotation angle can be adjusted by loosening nuts ³⁸ on the seat support flanges assembly ⁸.

13. For safety reasons, check the lower stud bolt ¹⁹. If necessary, it can be cut off in order to prevent rubbing.

Maintenance Tips

Diary

Rinse the lift with fresh pressurized water.

Perform a minimum of two up and down maneuvers to prevent deposits on the sleeve.

At the end of the day, close the water shut-off valve to release the pressure from the sleeve and avoid unnecessary stress on the seals.

Every week

Wipe the entire lift with a damp cloth to remove small lime deposits or residues.

Quarterly

Brush the supplied maintenance product over the entire steel structure of the elevator, let it act for a minimum of 3 hours (maximum 10 h), and rinse thoroughly with fresh pressurized water.

**The application of these tips is necessary to enable the warranty.*

**If superficial rust formations or many lime stains are observed, apply the supplied product for steel.*

Recommendations

During periods of inactivity exceeding 2 months, it is advisable to treat the steel with the supplied product (apply, let it act, and rinse). Once thoroughly dry, disassemble and store the elevator in a dry and protected place.

In the face of the danger of frost, it is advisable to protect and/or remove the elevator for shelter. For example, it can be wrapped with blankets or plastics.

When resuming the bathing season, ensure that the installation and fixation of the elevator are in perfect condition. Connect the power supply and perform 3-4 cycles to purge the air from the sleeve. (Follow steps 7 to 12 of the installation instructions).

The elevator should be annually inspected in its technical inspection by qualified personnel. A guide for this purpose will be provided. Fill out the technical inspection book located at the end of this manual.

Troubleshooting

WATER LEAKAGE

Possible causes:

Solutions:

Joint wear

Replace joints. See practical guide >>

Solid particles inside the cylinder

Disassemble the cylinder and clean all.
Install a filter in the network connection to prevent future damage.

Water overpressure

Install a reducer if the pressure is higher than 5 bars.

IT DOESN'T GO UP

Possible causes:

Solutions:

Low water supply pressure

Ensure that the pressure from the water network is within the regular working pressure range (3,5-5 bars).

Incorrect power supply

Make sure the power supply is properly connected and that the water passage is open.

Water inlet circuit clogged

Disassemble, clean all connectors and water inlets, both in the elevator and the external connector.

Incorrect handle position

Push the knob upwards until it reaches its stop position.

IT DOES NOT DESCEND

Possible causes:

Solutions:

Incorrect handle position

Push the knob upwards until it reaches its stop position.



Incrustations in the cylinder rail tube

Clean the translation rail, the shaft, and the cylinder. Install a filter in the network connection to prevent future damage.

Warnings

- > The hydraulic lift is a device strictly designed to facilitate water access for individuals with reduced mobility and must be used for this purpose.
- > The user must be alone and always in a seated position on the chair.
- > The lift should not be used with a weight exceeding the allowed limit (150 kg).
- > Do not force the operating lever beyond its travel.
- > Do not use the equipment with a network pressure exceeding 5 bars.
- > Use drinking water for its operation since the used water is returned to the pool during descent.
- > When not in use, always leave the chair out of the water. After each use, the elevator chair must always be raised.
- > Before using the elevator, make sure that the area of descent is clear.
- > Be cautious with the steel if the lift has been exposed to the sun for an extended period of time, as it may heat up.
- > Perform the necessary periodic maintenance and inspections to ensure the perfect operating conditions of the elevator.
- > Do not manipulate any part of the elevator unless by maintenance personnel or duly qualified for this purpose.

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